



Recombinant Schizosaccharomyces pombe Zinc finger protein sfp1 (sfp1)

Product Code	CSB-EP892418SXV-B
Abbreviation	sfp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UTH8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPSLALPINK PSHHNVNYNG NSFNSIHATS FGMSPQSWG N SFSGQAWLRD TIPLSNVVE SQTIP EEDST SYLNRLEEAF CRDFRCCGQT LEDLHQLIHH YEEQHAVLAT DSAVPQEYSL DSNNAAQSNH SHIQALKQRE RIRMHDLADQ LGTSEISDNS AVLPFAFPAN GGAPGPYRVS VVVPAAAAAA AAAASSDMSS DEASSQAETT GTPKKMPESL VMDASSPLSD MSMSIDVGES AANNVFAFNQ KDMVDSTYLP PFNYDHDVFS FAPSVASADQ FTESSMSPTP EVVSPAATNS AISSPFVRKS SSDLEAKPSK KQRSTPAFSH DSPLTIDYPG SNLVVVDKPY KCPVPNCDKA YKNQNGLYKH KLHGHCSPIT TPTPAPIPHQ GFVVENKPYR CEVCSKRYKN LNGLKYHRTH SHLQV SMAQA QREVQMNFM R TA
Source	E.coli
Target Names	sfp1
Protein Names	Recommended name: Zinc finger protein sfp1
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.