



# Recombinant Schizosaccharomyces pombe Zinc finger protein sfp1 (sfp1)

<b>Product Code</b>	CSB-MP892418SXV
<b>Abbreviation</b>	sfp1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UTH8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSLALPINK PSHHNVNYNG NSFNSIHATS FGMSPQSWG N SFSGQAWLRD TIPLSNVVE SQTIP EEDST SYLNRLEEAF CRDFRCCGQT LEDLHQLIHH YEEQHAVLAT DSAVPQEYSL DSNNAAQSNH SHIQALKQRE RIRMHDLADQ LGTSEISDNS AVLPFAFPAN GGAPGPYRVS VVVPAAAAAA AAAASSDMSS DEASSQAETT GTPKKMPESL VMDASSPLSD MSMSIDVGES AANNVFAFNQ KDMVDSTYLP PFNYDHDVFS FAPSVASADQ FTESSMSPTP EVVSPAATNS AISSPFVRKS SSDLEAKPSK KQRSTPAFSH DSPLTIDYPG SNLVVVDKPY KCPVPNCDKA YKNQNGLYKH KLHGHCSPIT TPTPAPIPHQ GFVVENKPYR CEVCSKRYKN LNGLKYHRTH SHLQV SMAQA QREVQMNFM R TA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	sfp1
<b>Protein Names</b>	Recommended name: Zinc finger protein sfp1
<b>Expression Region</b>	1-442
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.