



# Recombinant *Schizosaccharomyces pombe* Metaxin-like protein C409.19c (SPBC409.19c)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP892433SXV                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Abbreviation</b>      | SPBC409.19c                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | Q9UUA5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Immunogen Species</b> | <i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Sequence</b>          | MLQLFIYSPG LGQPTMDPGC LAALIYCALA VPKDEIEILR TANSGMSPTH<br>KLPALWDGHV WIGSLKNILI YLKQKGYNLD NFDKQLANV MAFTSLLEGS<br>VNDLWLEAF VNEENFVEAI RPAWSKALKF PHNYLTPNAL QRQAKERLAQ<br>TLGIRDEEVS YEASRMPISH KWTNATRHRQ ALLRTQARRI RISSLARQVY<br>GSLESLISDS KFIFGEKPTS LDCLFYAYLS FHAFTNELPQ ATRPCLQFN<br>SPKLYAYLKS LRETWFSDDS NILSPLSIKV QPENLLTIAR LAWNNVTAKA<br>NDTRKSITKF SVPPERKLLW ARNGFFIFAS AFSFVWFVIS NGIVVIETED<br>DEASFEEVDD AKESIQEKEI DETTESKATH DSSETSSSKE LPKEEEKESS<br>SFDLQPLSAQ DLLFSGFAED EIMDEEFGYD DDDDEEFDLD DLDDLEEEIV |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>      | SPBC409.19c                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Names</b>     | Recommended name: Metaxin-like protein C409.19c                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-450                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                         |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                  |