



Recombinant Human Cyclin-L1 (CCNL1)

Product Code	CSB-YP892449HU
Abbreviation	CCNL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UK58
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASGPHSTAT AAAAASSAAP SAGGSSSGTT TTTTTTTGGI LIGDRLYSEV SLTIDHSLIP EERLSPTPSM QDGLDLPSET DLRILGCELI QAAGILLRLP QVAMATGQVL FHRFFYSKSF VKHSFEIVAM ACINLASKIE EAPRRIRDVI NVFHHLRQLR GKRTPSPLIL DQNYINTKNQ VIKAERRVLK ELGFCVHVKH PHKIIVMYLQ VLECERNQTL VQTAWNYMND SLRTNVFVRF QPETIACACI YLAARALQIP LPTRPHWFLF FGTTEEEIQE ICIELRLTYT RKKPNYELLE KEVEKRKVAL QEAKLKAKGL NPDGTPALST LGGFSPASKP SSPREVKAAE KSPISINVKT VKKEPEDRQQ ASKSPYNGVR KDSKRSRNSR SASRSRSTR SRSRSHTPRR HYNRRSRSG TYSSRSRSRS RSHSESPRRH HNHGSPHLKA KHTRDDLKSS NRHGHKRKKS RSRSQSKSRD HSDAAKHRH ERGHHRDRRE RSRSFERSHK SKHHGGSRSRSG HGRHRR
Source	Yeast
Target Names	CCNL1
Protein Names	Recommended name: Cyclin-L1 Short name= Cyclin-L
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.