



Recombinant Human mRNA turnover protein 4 homolog (MRTO4)

Product Code	CSB-EP892452HU
Abbreviation	MRTO4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKD2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPKSKRDKKV SLTKTAKKGL ELKQNLIEEL RKCVDITYKYL FIFSVANMRN SKLKDIRNAW KHSRMFFGKN KVMMVALGRS PSDEYKDNLH QVSKRLRGEV GLLFTNRTKE EVNEWFTKYT EMDYARAGNK AAFTVSLDPG PLEQFPHSME PQLRQLGLPT ALKRGVVTLL SDYEVCKEGD VLTPEQARVL KLFQYEMAEF KVTIKYMWDS QSGRFQQMGD DLPESASEST EESDSEDDDD
Source	E.coli
Target Names	MRTO4
Protein Names	Recommended name: mRNA turnover protein 4 homolog
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein sharing a low level of sequence similarity with ribosomal protein P0. While the precise function of the encoded protein is currently unknown, it appears to be involved in mRNA turnover and ribosome assembly.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.