



Recombinant Human Protein kinase C and casein kinase substrate in neurons protein 3 (PACSIN3)

Product Code	CSB-MP892456HU
Abbreviation	PACSIN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKS6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPEEDAGGE ALGGSFWEAG NYRRTVQRVE DGHRLCGDLV SCFQERARIE KAYAQQQLADW ARKWRGTVEK GPQYGTLEKA WHAFFTAAER LSALHLEVRE KLQQQDSERV RAWQRGAFHR PVLGGFRESR AAEDGFRKAQ KPWLKRLKEV EASKKSYHAA RKDEKTAQTR ESHAKADSAV SQEQLRKLQE RVERCAKEAE KTKAQYEQTL AELHRYTPRY MEDMEQAFET CQAAERQRL FFKDMLLTLH QHLDLSSEK FHELHRDLHQ GIEAASDEED LRWWRSTHGP GMAMNWPQFE EWSLDTQRTI SRKEKGGRRSP DEVTLTSIVP TRDGTAPPPQ SPGSPGTGQD EEWSDEESPR KAATGVRVRA LYDYAGQEAD ELSFRAGEEL LKMSEEDEQG WCQGQLQSGR IGLYPANYVE CVGA
Source	Mammalian cell
Target Names	PACSIN3
Protein Names	Recommended name: Protein kinase C and casein kinase substrate in neurons protein 3 Alternative name(s): Endophilin I SH3 domain-containing protein 6511
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.