



# Recombinant Human Protein kinase C and casein kinase substrate in neurons protein 3 (PACSIN3)

<b>Product Code</b>	CSB-YP892456HU
<b>Abbreviation</b>	PACSIN3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UKS6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAPEEDAGGE ALGGSFWEAG NYRRTVQRVE DGHRLCGDLV SCFQERARIE KAYAQQQLADW ARKWRGTVEK GPQYGTLEKA WHAFFTAAER LSALHLEVRE KLQQQDSERV RAWQRGAFHR PVLGGFRESR AAEDGFRKAQ KPWLKRLKEV EASKKSYHAA RKDEKTAQTR ESHAKADSAV SQEQLRKLQE RVERCAKEAE KTKAQYEQTL AELHRYTPRY MEDMEQAFET CQAAERQRL FFKDMLLTLH QHLDLSSEK FHELHRDLHQ GIEAASDEED LRWWRSTHGP GMAMNWPQFE EWSLDTQRTI SRKEKGGRRSP DEVTLTSIVP TRDGTAPPPQ SPGSPGTGQD EEWSDEESPR KAATGVRVRA LYDYAGQEAD ELSFRAGEEL LKMSEEDEQG WCQGQLQSGR IGLYPANYVE CVGA
<b>Source</b>	Yeast
<b>Target Names</b>	PACSIN3
<b>Protein Names</b>	Recommended name: Protein kinase C and casein kinase substrate in neurons protein 3 Alternative name(s): Endophilin I SH3 domain-containing protein 6511
<b>Expression Region</b>	1-424
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.