



# Recombinant Human Isobutyryl-CoA dehydrogenase, mitochondrial (ACAD8)

<b>Product Code</b>	CSB-EP892457HU-B
<b>Abbreviation</b>	ACAD8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UKU7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LVQTGHRSLTSCIDPSMGLNEEQKEFQK VAFDFAAREM APNMAEWDQK ELFPVDVMRK AAQLGFGGVY IQTDVGGSGLSRLDTSVIFE ALATGCTSTT AYISIHNMCA WMIDSFGNEE QRHKFCPPLC TMEKFASYCL TEPGSGSDAA SLLTSAKKQG DHYILNGSKA FISGAGESDI YVVMCRTGGP GPKGISCIVV EKGTPGLSFG KKEKKVGWNS QPTRAVIFED CAVPVANRIG SEGQGFLIAV RGLNGGRINI ASCSLGAAHA SVILTRDHLN VRKQFGEPLA SNQYLQFTLA DMATRLVAAR LMVRNAVAL QEERKDAVAL CSMAKLFATD ECFAICNQAL QMHGGYGYLK DYAVQQYVRD SRVHQILEGS NEVMRILISR SLLQE
<b>Source</b>	E.coli
<b>Target Names</b>	ACAD8
<b>Protein Names</b>	Recommended name: Isobutyryl-CoA dehydrogenase, mitochondrial EC= 1.3.99.- Alternative name(s): Activator-recruited cofactor 42 kDa component Short name= ARC42 Acyl-CoA dehydrogenase family member 8 Short name= ACAD-8
<b>Expression Region</b>	23-415
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.