



Recombinant Human Coiled-coil domain-containing protein 19, mitochondrial (CCDC19)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-BP892460HU |
| Abbreviation | CCDC19 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9UL16 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | KAV SSEVDESLFG DIKSPAQQQS DSPIVLLRDK HTLQKTLTAL GLDRKPETIQ LITRDMVREL IVPTEDPSGE SLIISPEEFE RIKWASHVLT REELEARDQA FKKEKEATMD AVMTRKKIMK QKEMVWNNNK KLSDLEEVAK ERAQNLLQRA NKLRMEQEEE LKDMSKIILN AKCHAIRDAQ ILEKQQIQKE LDTEEKRLDQ MMEVERQKSI QRQEELERKR REERIRGRRQ IVEQMEKNQE ERSLAEQRE QEKEQMLEYM EQLQEEDLKD MERRQQQKLK MQAEIKRIND ENQKQKAELL AQEKLADQMV MEFTKKKMAR EAEFEAEQER IRREKEKEIA RLRAMQEKAQ DYQAEQDALR AKRNQEVADR EWRRKEKENA RKKMETEAEL RKSRLQVAF KEHALAVQVQ RDRDEFERIL RAQREQIEKE RLEEEKKATG RLQHANELRR QVRENQQKEV QNRIATFEEG RRLKEEAQKR RERIDEIKRK KLEELRATGL PEKYCIEAER KANILPATSV N |
| Source | Baculovirus |
| Target Names | CFAP45 |
| Protein Names | Recommended name: Coiled-coil domain-containing protein 19, mitochondrial Alternative name(s): Nasopharyngeal epithelium-specific protein 1 |
| Expression Region | 28-551 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.