



Recombinant Human Coiled-coil domain-containing protein 19, mitochondrial (CCDC19)

Product Code	CSB-YP892460HU
Abbreviation	CCDC19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UL16
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	KAV SSEVDESFLG DIKSPAQQQS DSPIVLLRDK HTLQKTLTAL GLDRKPETIQ LITRDMVREL IVPTEDPSGE SLIISPEEFE RIKWASHVLT REELEARDQA FKKEKEATMD AVMTRKKIMK QKEMVWNNNK KLSDLEEVAK ERAQNLLQRA NKLRMEQEEE LKDMSKIILN AKCHAIRDAQ ILEKQQIQKE LDTEEKRLDQ MMEVERQKSI QRQEELERKR REERIRGRRQ IVEQMEKNQE ERSLLAEQRE QEKEQMLEYM EQLQEEDLKD MERRQQQKLK MQAEIKRIND ENQKQKAELL AQEKLADQMV MEFTKKKMAR EAEFEAEQER IRREKEKEIA RLRAMQEKAQ DYQAEQDALR AKRNQEVADR EWRRKEKENA RKKMETEAEEL RKSRLQVAF KEHALAVQVQ RDRDEFERIL RAQREQIEKE RLEEEKKATG RLQHANELRR QVRENQQKEV QNRIATFEEG RRLKEEAQKR RERIDEIKRK KLEELRATGL PEKYCIEAER KANILPATSV N
Source	Yeast
Target Names	CFAP45
Protein Names	Recommended name: Coiled-coil domain-containing protein 19, mitochondrial Alternative name(s): Nasopharyngeal epithelium-specific protein 1
Expression Region	28-551
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.