



Recombinant Human RNA-binding protein NOB1 (NOB1)

Product Code	CSB-EP892472HU-B
Abbreviation	NOB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULX3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	APVEHVVD AGAFLRHAAL QDIGKNIYTI REVVTEIRDK ATRRRLAVLP YELRFKEPLP EYVRLVTEFS KKTGDYPSLS ATDIQVLALT YQLEAEFVGV SHLKQEPQKV KVSSSIQHPE TPLHISGFHL PYKPKPPQET EKGHSACEPE NLEFSSFMFV RNPLPNIDHE LQELLIDRGE DVPSEEEEEEE ENGFEDRKDD SDDGGGGWIT PSNIKQIQQE LEQCDVPEDV RVGCLTTDFA MQNVLLQMGL HVLAVNGMLI REARSYILRC HGCFKTTSDM SRVFCSHCGN KTLKKVSVTV SDDGTLHMHF SRNPKVLNPR GLRYSLPTPK GGYAINPHL TEDQRFQQLR LSQKARQKTN VFAPDYIAGV SPFVENDISS RSATLQVRDS TLGAGRRRLN PNASRKKFVK KR
Source	E.coli
Target Names	NOB1
Protein Names	Recommended name: RNA-binding protein NOB1 Alternative name(s): Phosphorylation regulatory protein HP-10 Protein ART-4
Expression Region	2-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.