



Recombinant Human Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide (DAPP1)

Product Code	CSB-BP892478HU
Abbreviation	DAPP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UN19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGRAELLEGK MSTQDPSDLW SRSDGEAELL QDLGWYHGNI TRHAAEALL SNGCDGSYLL RDSNETTGLY SLSVRAKDSV KHFHVEYTG YSFKFGFNEFS SLKDFVKHFA NQPLIGSETG TLMVLKHPYP RKVEEPSIYE SVRVHTAMQT GRTEDDLVP TAPSLGTKEGY LTKQGGLVKT WKTRWFTLHR NELKYFKDQM SPEPIRILD LTECSAVQFDY SQERVNCFCL VFPFRTFYLC AKTGVEADEW IKILRWKLSQ IRKQLNQGEG TIRSRSFIFK
Source	Baculovirus
Target Names	DAPP1
Protein Names	Recommended name: Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide Short name= hDAPP1 Alternative name(s): B lymphocyte adapter protein Bam32 B-cell adapter molecule of 32 kDa
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.