



Recombinant Schizosaccharomyces pombe Cell wall synthesis protein psu1 (psu1)

Product Code	CSB-MP892500SXV
Abbreviation	psu1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UR09
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	AP FRHPHLLNK RDVSVVTSKV YAYTTVTLEA AASAISTNGA AKEAATAAGD AETTSSVAAS VTPAASSSVA ASVTPVASSS VAASVTPVSS SAVDSATSA AASSSVIPTS SSVVASSEV ASSTTSSAAA SATSTGSSSG GFQDGVYDCT DFPDQSGVV ALDYLGYYGGY SGIQIGDGAG SSCVEGAYCS YACSPGMLKT QWPSTQPSDG ETRGLLCKG GKLYRTNTAY DNLCEGVGT ASVQNTLGQG VAICQTDYPG SENMVIPTYV GGGATSPLSV TDNANYFQWE GKG TSAQYYV NKAGYTAEQG CQWGTSSGDY GNWSPLVFGA GMTDGTWLS ISQNPLTSQL ANYNVKIVGA DGASVSGTCV FENGAFQNGG SGCTVGITSG SGVVFVY
Source	Mammalian cell
Target Names	psu1
Protein Names	Recommended name: Cell wall synthesis protein psu1
Expression Region	19-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.