



Recombinant *Saccharomyces cerevisiae* tRNA-specific adenosine deaminase subunit TAD3 (TAD3)

Product Code	CSB-EP892502SVG-B
Abbreviation	TAD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9URQ3
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVKKVNNPLK IDYQNGIIE N RLLQIRNFKD VNTPKLINVW SIRIDPRDSK KVIELIRNDF QKNDPVSLRH LKRIRKDIET STLEVVLCSK EYICDEGEIN NKLKSIWVGT KKYELSDDIE VPEFAPSTKE LNNAWSVKYW PLIWNNGNPND QILNDYKIDM QEVNLSRA STLSVKMATA GKQFPMVSF VDP SRKKDKV VAEDGRNCEN SLPIDHSVMV GIRAVGERLR EGVDEDANSY LCLDYDVYLT HEPCSMCSMA LIHSRVRV V FLTEMQRTGS LKLTSGDGYC MNDNKQLNST YEAFAQWIGEE YPVGQVDRDV CC
Source	E.coli
Target Names	TAD3
Protein Names	Recommended name: tRNA-specific adenosine deaminase subunit TAD3 Alternative name(s): tRNA-specific adenosine-34 deaminase subunit TAD3
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.