



Recombinant Schizosaccharomyces pombe Uncharacterized protein C4B3.18 (SPCC4B3.18)

Product Code	CSB-BP892522SXV
Abbreviation	SPCC4B3.18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USK7
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPEMDKSNDF FEENPCPNL GDVSSLIEDF INLQSNQNI AFVTSGGTLV PLEQNTVRFI DNFSAGNRGA ASAEYFCKSG YSVIFLYRNY SLMPFTRRYE GPWTNLFERS NGSQEKPWIV KQEFQSTVFN ALNELDTYTK HHRLMIPYT TLTEYLWYLK AIAERLQIIA SRALFYLA AAA VSDFHIPPKE LSEHKIQSGS GSNKLEVKMH PVPKFLRFLV DTWAPEASII SFKLETDP SI LIDKARVALQ RYNHQLVIAN LLSTRKRSVA FVTSKTVEWL HLSDDDDIVKK VEIEQYIVSR AISLHKERIQ EFAKAP
Source	Baculovirus
Target Names	SPCC4B3.18
Protein Names	Recommended name: Uncharacterized protein C4B3.18
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.