



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C1805.14 (SPAC1805.14)

Product Code	CSB-BP892553SXV
Abbreviation	SPAC1805.14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UTG3
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFRRNTLTPG KGRNSLDMFS VNDFMTWIDG MKKTSDENGC SPFTLKPIKS TSDRKRRRHS RMSIVNVWNQ SPTSYRDQLL HFNHESDVSI DDTSEVAGEE TFGQIERFEY LLSSDGDGDG DVEIENPEVD DYQVSEHDFG LEDDFNGDGC QQMIEIESDE AISEEEEDLS NENKSESQGL ESFSFNQDNT FVNYASPASI TEEFVEPKDS QFNHDVNLQM GSAPGYSTVE PEDNFASEIQ TNAPEVHLNY ENSDYTEDI DLLDHHFCDL SEISKFNHQH SGKPDHPSLV SNASLAPFIV EGNGIKNGLL HYNMETAETD ESYTDLDDTF ARLKNLSQRH TNHSTDHDD TVDSHLLHSF VSAENANEPT NDVDDNSLQE QVADASQFVS FLETTKTIVAT SNLRSTKRKL SEILENGQSG DCSLTDTNVD FDISEKQASG NSRSMIPLRK KR
Source	Baculovirus
Target Names	SPAC1805.14
Protein Names	Recommended name: Uncharacterized protein C1805.14
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.