



# Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C1805.14 (SPAC1805.14)

<b>Product Code</b>	CSB-MP892553SXV
<b>Abbreviation</b>	SPAC1805.14
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UTG3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFRRNTLTPG KGRNSLDMFS VNDFMTWIDG MKKTSDENGC SPFTLKPIKS TSDRKRRRHS RMSIVNVWNQ SPTSRYRDQLL HFNHESDVSI DDTSEVAGEE TFGQIERFEY LLSSDGDGDG DVEIENPEVD DYQVSEHDFG LEDDFNGDGC QQMIEIESDE AISEEEEDLS NENKSESQGL ESFSFNQDNT FVNYASPASI TEEFVEPKDS QFNHDVNLQM GSAPGYSTVE PEDNFASEIQ TNAPEVHLNY ENSDYTEDHI DLLDHHFCDL SEISKFNHQH SGKPDHPSLV SNASLAPFIV EGNGIKNGLL HYNMETAETD ESYTDLDDTF ARLKNLSQRH TNHSTDHDD TVDSHLLHSF VSAENANEPT NDVDDNSLQE QVADASQFVS FLETTKT VAT SNLRSTKRKL SEILENGQSG DCSLTDTNVD FDISEKQASG NSRSMIPLRK KR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPAC1805.14
<b>Protein Names</b>	Recommended name: Uncharacterized protein C1805.14
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.