



Recombinant Schizosaccharomyces pombe Serine/threonine-protein kinase ppk34 (ppk34)

Product Code	CSB-MP892568SXV
Abbreviation	ppk34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UU87
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSLPLPSFLS SNNYNSSETS LLTEGPESDE EDEGPLLKQY RLKNMLGYGA CSTVYLAVDV STNIEYAIKE FKKTSLRRRE KFRLMQLEKE LGSFNDMDSS TIEDQIRNQE KKDPYYYIRK ELDVLKVLDH ENVVKLYQVI NSDYKDSL VLNYCPGGKL ASVDHGISFD PMPLETCKRS FRQLVLAVNY IHGKGVHRD IKPDNLFKE NNNLCVIDFG LAEFLPKDGF VKPASGTPAF MAPELFDNEL KKIKGKPLDI WSMGVTLTYVI AFAEFPFNGS TILDMINVIK GKSLMIPAYC NSDLRKLLEK CLEKNPEKRI TIEELLRDPF LTERGKHHLT SSSIVTAFSS IWEI
Source	Mammalian cell
Target Names	ppk34
Protein Names	Recommended name: Serine/threonine-protein kinase ppk34 EC= 2.7.11.1
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.