



Recombinant *Pyrococcus abyssi* Cell division control protein 6 homolog (cdc6)

Product Code	CSB-MP892667FHV
Abbreviation	cdc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V2F2
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrococcus abyssi</i> (strain GE5 / Orsay)
Purity	≥85% (SDS-PAGE)
Sequence	MSEGEQLHLD KLFERLLKAR KIFKNKEVLR HSYTPKDLPH RHEQIEALAQ ILVPVLKGET PSNIFVYGKT GTGKTVTVKF VTEELKKVSQ KYNIPVEVIY INCEIVDTHY RVLANIVNHF KDETGIEVPL VGWPTDEVYA KLKQVIDMKE RFVIIVLDEI DKLVKKSGDE VLYSLTRINT ELKRAKVSVI GISNDLKFKKE YLDPRVLSSL SEEEVVFPPY DANQLKDILT QRAEEAFYPG VLEEGVIPLC AALAAREHGD ARKALDLLRV AGEIAEREGA SKVTEKHVWR AQEKIEQDMM EEVIKTLPLQ SKVLLYAILV LDENGDLPAN TGEVYSVYRE LCEYLDLEPL TQRRISDLIN ELDMLGIINA KVVSCKGRYGR TKEIRLNVTP YKIRNVFRYD YAIQPLLTVS LKGGQRLL
Source	Mammalian cell
Target Names	cdc6
Protein Names	Recommended name: Cell division control protein 6 homolog Short name= CDC6 homolog
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.