



Recombinant *Drosophila melanogaster* Ero1-like protein (Ero1L)

Product Code	CSB-MP892675DLU
Abbreviation	Ero1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V3A6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	G YFAAIDETET SKNCFCELEG SINDCSCDVD TVDHFNNMKI YPRLQSLLVK NFFRFYKVNLRQECFPWPDD SRCAMRFCQV ENCEEQAIPQ GIKDKGEHKE KAAFKYTREAA QVGGSSACSDG EDFDSSLGFL DTSISDQAGR EFELWAKHDE AEEDFCIVDD HEEGSQYVDL LLNPERYTG Y KGESAHRIWK SIYLENCFGG NNETANKFSN YVPHLDLRNV CLEQRAF YRI ISGLHSSINI HLCSKYLLSE SKDFLDPQGI WGPNVKEFKR RFSPETTSGE GPHWLRNLYF IYLIELRALA KAAPYLRRED YYTGIAEEDD EVKLAIN DML SVIENFQSHF DENALFSNGI ASIKFKHDYK EKFRNISRIM SCVGC DKCKL WGKLQTQGLG TALKILYSEK LNLATESGLW DKPHIEADPI FRLSRTEIVA LFNAFGRLSN SIYEMENFRC VLR
Source	Mammalian cell
Target Names	Ero1L
Protein Names	Recommended name: Ero1-like protein EC= 1.8.4.- Alternative name(s): Endoplasmic oxidoreductin-1-like protein
Expression Region	30-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.