



Recombinant *Drosophila melanogaster* Kinetochores protein Spc25 (Spc25)

Product Code	CSB-BP892679DLU
Abbreviation	Spc25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V3V7
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAIIMTESSY ERRVKALY EK QIHMEALEAK FIKKVFKEFNS NLLDVKEAAS RHQRKVGK LQ KVLMERREEL DKRVSFIEEL DRELEATKLH NLAMKDWFKQ QKMLAKQRKN EIMESIHTLS KTRTYINQE ALPARVKGVT VLRGDKRDQL IPFDLKATDV EGLDSLQHL ESLNVDVAQW QQLISLAMDM AMESRAPPTTP PKEADNCKSI IEIDLTSPMS HT
Source	Baculovirus
Target Names	Spc25
Protein Names	Recommended name: Kinetochores protein Spc25 Alternative name(s): Chromosome alignment defect 2 Mitotic chaos
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.