



Recombinant *Drosophila melanogaster* Putative proteasome inhibitor (CG8979)

Product Code	CSB-EP892687DLU
Abbreviation	CG8979
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V637
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MESGSTAKTG DFFYGWDL LY KTVKADVSKK SDLLIALVHF LLTKHYNFRC VGVGDDKTLP EEESSELLPD SWNDDDTKYS LRYVHDKMLY LLLGHITEGS LLINLLDINT KKVSNICVEP ETLVPEVKGG ITTIMPSASE IVERYRRELL DPVFTGNSRE VTTQTTNSPR PIGSDPDPLR IGEPRRGGSF IPSAFEPRPF GFPDVGRGDL DPLGRGGHGN LFSFSPRPNM GPGPVPRFDP FNPLNPNRPG QGGINPDHMR PPDWNPDYIM
Source	E.coli
Target Names	PI31
Protein Names	Recommended name: Putative proteasome inhibitor
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.