



Recombinant *Drosophila melanogaster* Acidic leucine-rich nuclear phosphoprotein 32 family member A (Anp32a)

Product Code	CSB-EP892694DLU-B
Abbreviation	Anp32a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V895
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEKRIELERR ARKVNQITEL NLDNCRSTSI VGLTDEYTALESLSLINVGL TTLKGFPKLP NLKKLELSDN RISSGLNYLT TSPKLQYLNL SGNKIKDLET LKPLEEFKNL VVLDLFNDA TQVDNYREKI FKMLPSLNFL DGFDCNDEEV QSDGDDDDDEV NGNDSDEVGVSDEDDDSDDS DEEANGEVSL SEVYNDLEE DNSDWEGEDE AGEEDDEEEDS DIDDADGDAN ESAASVNAKD KDGEKEADES QVRGKKRKHD G
Source	E.coli
Target Names	Anp32a
Protein Names	Recommended name: Acidic leucine-rich nuclear phosphoprotein 32 family member A
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.