



Recombinant *Drosophila melanogaster* Spindle assembly abnormal protein 6 homolog (sas-6)

Product Code	CSB-BP892701DLU
Abbreviation	sas-6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VAC8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MWPPGSEDSY SAKMDYGKSV VNILPSVEML VNFNGEMTRS SKRSCLLYAE RVDFKELLQL RLTEKSDQRR MYITTVDSAS FQDLKQDQSL NVSFSGFIDN VVRMLKDCQS GKLELHLTTR DQNLSSGREV HDYYLQFVEI RSFKNLVHLS LPCRSAPLNT VLFYINSMLE ASHKKQYILE QSMQQMQAEI NAQRAHAERL TTENTNIREA LAENTRILEE KHAAEVHQYQ EKLSKINEQR SNELENNRRA ISGFQAQLDK ASLEKAEKLS AQEQAERKRCQ TLSEELSCCK ARVCTLKEQN DKLHG DVANI RKHERKLEYK IEDLKQHTVE LQEHQKGNK EKANIAAELE AEKKILHTKR QALEMASEEI SKANQIIVKQ SQELLNLKKT IAWRTEVALQ QEKAVQAKES LLSLRENELR EARITIEKLR EEIPQQLQSM RNFAQGLEQK YSKQILILKE RLAIPTGKEN RR
Source	Baculovirus
Target Names	Sas-6
Protein Names	Recommended name: Spindle assembly abnormal protein 6 homolog
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.