



Recombinant Drosophila melanogaster Zinc finger protein 330 homolog (Noa36)

Product Code	CSB-YP892703DLU
Abbreviation	Noa36
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VAU9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPKKKTGQRK KAEKQKLRLK EIRSREVPLA DLPCNAPMEC DKCEKKQKSR AFCYFCQSIQ RLPICAQCGK IKCMLKTGDC VVKHPGVYTT GLGMVGAICD FCEAWVCHGR KCLQNHACTC PLQNATCLEC ERGVWEHGGR IFKCSFCNGF LCEDDQFEHQ ASCQVLESEN YKCQSCNKLK QYSCLRCKTC YCEDHVRKRG FKYDKNKPIP CPKCNYDTSV TKDLSMSTRS HKFGRQQQGG NSDDEEGYGG YYGAGGSGYY GGAASGGYSY GGDDDEDESD GDYDDESDDED DDDDEEEDDET TESEPEKETD KKATK
Source	Yeast
Target Names	Noa36
Protein Names	Recommended name: Zinc finger protein 330 homolog Alternative name(s): Nucleolar autoantigen 36 homolog
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.