



Recombinant Drosophila melanogaster Glutathione S-transferase D2 (GstD2)

Product Code	CSB-BP892842DLU
Abbreviation	GstD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VG98
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDFYYMPGGG GCRTVIMVAK ALGLELNKKL LNTMEGEQLK PEFVKLN PQH TIPTLVDNGF SIWESRAIAV YLVEKYGKDD YLLPNDPKKR AVINQRLYFD MGTLYESFAK YYYPLFRTGK PGSEDEDLKRI ETAFGLDTF LEGQEYVAGD QLTVADIAIL STVSTFEVSE FDFSKYSNVS RWYDNAKKVT PGWDENWEGL MAMKALFDAR KLA AK
Source	Baculovirus
Target Names	GstD2
Protein Names	Recommended name: Glutathione S-transferase D2 Short name= DmGST21 EC= 2.5.1.18
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.