



Recombinant Drosophila melanogaster Protein

Fer3 (fer3)

Product Code	CSB-EP892844DLU-B
Abbreviation	fer3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VGJ5
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MQHPHPIDQP TYMPDVFPFQP LWGQEAPPPP IVPYQELIAG FPCTDLSLWQ RSQVTPLVPQ RPSTNGRANG SSSSSKTRR RVASMAQRRR ANIRERRRMF NLNEAFDKLR RKVPTFAYEK RLSRIETLRL AITYIGFMAE LLSGTPSNSH KSRSDVYGSM NGHHQAPPPA IHPHHLHPAA AYQRDFASPY NHSLS
Source	E.coli
Target Names	fer3
Protein Names	Recommended name: Protein Fer3 Alternative name(s): Basic helix-loop-helix protein N-twist Nephew of atonal 3 Neuronal twist
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.