



Recombinant Drosophila melanogaster Beta-lactamase-like protein 2 homolog (CG12375)

Product Code	CSB-EP892866DLU
Abbreviation	CG12375
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VLS9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MALIPAVTRL TSSVIRILGC NPSAMTLQGT NTYLLGSGSR RILIDTGDED VPQYIAHLGD VLQQEKASID TILLTHWHHD HVGGVKSIVG TKLAEKDCRV FKFGRTDAPD VCPEIPTDIK LHPLAHNQEFG TTEGANVRV VHTPGHTTDHV VLAMNEGTLF SGDCILGEGT AVFEDLFEYM KSLEKILDIK PQRIFPGHGN VIEEPIGKIE YYINHRNQRE QQILQFFVQR PNENLQAMDV VKVVYKETPE NLWPAAAYNV NHHLKLEKE GKLRVSKSAD EEFYVYQPTS AL
Source	E.coli
Target Names	CG12375
Protein Names	Recommended name: Beta-lactamase-like protein 2 homolog EC= 3.-.-.-
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.