



Recombinant *Drosophila melanogaster* RNA-binding protein pno1 (I (1)G0004)

Product Code	CSB-EP892885DLU-B
Abbreviation	I(1)G0004
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VR89
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEAENIRADA FEPAKRATKR GASGGGQQDV EMQVDEATGI EGQVLGSSRA SAPPKAKRAR SELRKVSVPP HRYSSLKEHW MKIFTPVEH MKLQIRFNMK ARQVELRVGP ETPDIANLQR GADFVRAFLC GFEVDDALAL LRLEDLFVES FEIKDVKTLR GDHQSRAIGR LAGKGGRTKF TIENVTKTRI VLADSKIHIL GSYQNIQLAR RAVCNLILGS PPSKVYGNLR AVASRLSERM
Source	E.coli
Target Names	I(1)G0004
Protein Names	Recommended name: RNA-binding protein pno1
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.