



Recombinant *Drosophila melanogaster* Probable insulin-like peptide 3 (Ilp3)

Product Code	CSB-YP892894DLU
Abbreviation	Ilp3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VT52
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	T MKLCGRKLPE TSKLCVYGF NAMTKRTLDP VNFNQIDGFE DRSLLELLS DSSVQMLKTR RLRDGVFDEC CLKSCTMDEV LRYCAAKPRT
Source	Yeast
Target Names	Ilp3
Protein Names	Recommended name: Probable insulin-like peptide 3 Short name= dILP3 Alternative name(s): Insulin-related peptide 3 Cleaved into the following 2 chains: 1. Probable insulin-like peptide 3 A chain 2. Probable insulin-like peptide 3 B
Expression Region	30-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.