



Recombinant Drosophila melanogaster PCI domain-containing protein 2 homolog (PCID2)

Product Code	CSB-MP892896DLU
Abbreviation	PCID2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VTL1
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MFGTVNNYLS GVLHAAQDLD GESLATYLSL RDVHVQNHNL YIAQPEKLVD RFLKPPLDEV VSAHLKVLHY LAQEPPGYME AYTQQSAACG AVVRLQQLK DENWCLPLMY RVCLDLRYLA QACEKHCQGF TPGHVLEKAA DCIMACFRVC AADGRASEED TKRLGMMNLV NQLFKIYFRI NKLHLCKPLI RAIDNCIFKD SFPLPEQITY KYFVGRRAMF DSNYQAAVQY LSYAFSNCPD RFASNKRLIL IYLVPVKMLL GYLPSKSLQ RYDLLLFLDL AMAMKAGNVN RFDEIVRDQE LVLIRSGIYL LVEKLFVY RNLFKKVFI RKSHQLDMGD FLSALHFVGL TDVSLDETHC IVANLIYDGK IKGYISHAHN KLVVSKQNPFSVSL
Source	Mammalian cell
Target Names	PCID2
Protein Names	Recommended name: PCI domain-containing protein 2 homolog Short name=dmPCID2 Alternative name(s): CSN12-like protein
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.