



Recombinant *Drosophila melanogaster* Transcription initiation factor TFIID subunit 8 (Taf8)

Product Code	CSB-EP892911DLU-B
Abbreviation	Taf8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VWY6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEKVEAVTVS NVDPYRRILN KVVSQLLLDK GAGQASNHSL ETLTQMLQAL IWEIGNSAHN YCELSGRTMP TVGDVSLALI NMGISISNLD PYMRKETHVP IPLPPQQTQQ RPLSLLQAGI KAPHPHYVPS YFPPMPDPHA YIRTPTHKQP VTEYEAIKREK AACQKRDIK ALTKFLCKTT ETNNLFPTED NMFPLIACKP AFPPYAAALN PTDQVDFDFFEE LEYHYLVANR TEDEPSKDDG EEGDSENEEM DGDKSKEEKP ELDIKPNSNT NKAILNPNI DNPYLRAATL PKRSKNCPPT GTMPRSRLAT TAPTIRTPST LEITKTNL
Source	<i>E.coli</i>
Target Names	Taf8
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 8 Alternative name(s): Protein prodos
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.