



Recombinant *Drosophila melanogaster* Uncharacterized protein C21orf59 homolog (CG18675)

Product Code	CSB-MP892921DLU
Abbreviation	CG18675
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VZH1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVLLQVKRGD ELIFLYETSV KEKTDTVLRE LVALHNGQLK IQRVCMEIEE LAEHGTMVPP EMIGLNEEQR EELKLKDVWA DKCIPSGGFT INKDLLRRN GQQPTEAMQK VLAGAMNDAK AMVDRRLAKS SKTLTPKIVE EALNLLRGGV TIVYPMQLPP HDSIRMEFTN TEDLTGTQAS KEVIEPSKAQ FWFAGHQMLM GKLMSDYLGH NDKTKVVVKL NQIGEGPPGR EAVISDPLRR QMMAEAYRRQ EELKKLEQDD DDEYLNSSWA DSGSLKRQAH GLENVRFRRFR K
Source	Mammalian cell
Target Names	CG18675
Protein Names	Recommended name: Uncharacterized protein C21orf59 homolog
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.