



Recombinant *Drosophila melanogaster* Chitinase-like protein Idgf4 (Idgf4)

Product Code	CSB-BP892933DLU
Abbreviation	Idgf4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W303
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	AGSHHLLCY YDGNFVREG LSKLILTDLE PALQYCTHLV YGYAGINPSS NKLVSNNKL DLDLGSSLFR QVTGLKRKYP ALKVLLSVGG DKDTPDPENN KYLTLLESSN ARIPFNSAH SLVKTYGFDG LDLGWQFPKN KPKKVHGSIG KFWKGFKKIF SGDHHVDEKA EEHKEAFTAL VRELKNAFRP DGYILGLSVL PNVNSSLFFD VPAIINLDY VNLHTYDFQT PERNNEVADF PAPIYELNER NPEFNVNYQV KYWTGNRAPA AKINVGIATY GRAWKLTKDS GLTGLPPVAE ADGVAPAGTQ TQIPGLLSWP EVCAKLPNPA NQHLKGADGP LRKVGDPDKR FGSYAYRSAD DSGENGVWVG YEDPDAAIK AEYVKREGLG GIAVVDLSFD DFRGGCTGHD KFPILRQVKS KL
Source	Baculovirus
Target Names	Idgf4
Protein Names	Recommended name: Chitinase-like protein Idgf4 Alternative name(s): Imaginal disk growth factor protein 4
Expression Region	22-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.