



Recombinant *Drosophila melanogaster* tRNA (guanine-N (7)-)-methyltransferase subunit wuho (wuho)

Product Code	CSB-EP892937DLU
Abbreviation	wuho
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W415
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MCTTISFAEP EIVLGHGRRV LFNVPDDLQI FKEIELPPDL GLKGHTSQSQ ESCTAAAAAS TATAASGQAP GGKEQQLANQ PEEGGTSASA SGLGCATSTS VQNVAYSPDG QLLAVTTSGK QKALLLYRSR PENARLLSAR PLARASSALR FCSDSSSILV TDKTGDCYQY DCVEVEAPPR LLLGHLSVVY DILWSEDQQH IITCDRDDKI RVTNYPATFD IHSYCLGHRE FVSGLLALTE QHIASASGDK TLRVWNYIQG KELLQHELPA PAVRLLVRQL EPEKVFQAAV LFYEHVDALG LYRLERSDD TWSVTATQLV CAEAGSWSIS NFTLTSDRIY ITGAENERLS LRVYDIATGQ PASSGVPEGW LKMVLDGLGA NEEGAPPFIP EDLSVWFKKR FDNVSDYLER KKRRIEEQQQ QKCG
Source	E.coli
Target Names	wuho
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase subunit wuho
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.