



Recombinant *Drosophila melanogaster* V-type proton ATPase subunit d 1 (VhaAC39-1)

Product Code	CSB-MP892938DLU
Abbreviation	VhaAC39-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W4P5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MNSSGFMFNI DNGYLEGLCR GFKCGILKQA DYLNLVQCET LEDLKLHLQG TDYGSFLANE PSPLSVSVID DKLREKLVIE FQHMRNHAVE PLSNFLDFIT YGYMIDNIIL LITGTLHQRP ISELIPKCHP LGSFEQMEAI HVASTPAELY NAVLVDTPLA PFFVDCISEQ DLDEMNIIEI RNTLYKAYLE AFYNFCKNMG GATADVMCEI LAFEADRAI IITINSFGTE LSKDDRACKLY PNCGKMYPDG LAALARADDY EQVKTVAEYY AEYAALFDGS GNNPGDKTLE DKFFEHEVKL DVYAFLQQFH FGVFYAYLKL KEQEARNIVW IAECVAQKHR AKIDNYIPIF
Source	Mammalian cell
Target Names	VhaAC39-1
Protein Names	Recommended name: V-type proton ATPase subunit d 1 Short name= V-ATPase subunit d 1 Alternative name(s): V-ATPase 39 kDa subunit Vacuolar H ⁺ ATPase subunit AC39-1 Vacuolar proton pump subunit d 1
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.