



Recombinant Chicken Proheparin-binding EGF-like growth factor (HBEGF), partial

Product Code	CSB-EP892955CH1
Abbreviation	HBEGF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W7C5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	AV GKFGRDGPFG EALHNEVLLH EGGGVASVSA TAPLLGGIVE KEGGGAASGD ALSELPRVAF LSKPQGPVTP KKKGNGNKRR KGKGLGKKRD PCLRKYKDFC IHGECKYIRE LGAPSCICQP GYHGERCHGL LLPVEHPPST YDHTTAL
Source	E.coli
Target Names	HBEGF
Protein Names	Recommended name: Proheparin-binding EGF-like growth factor Cleaved into the following chain: 1. Heparin-binding EGF-like growth factor Short name= 2. HB-EGF Short name= 3. HBEGF
Expression Region	19-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.