



Recombinant Mouse Origin recognition complex subunit 5 (Orc5)

Product Code	CSB-EP892993MO-B
Abbreviation	Orc5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUV0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSHLESMVLC REAQVSTLQS LFGERHHFSF PSIFIYGHTA SGKTYVTQTL LKTLELPHAF VNCVECFTSR LLEQILNKL SHLNSSDAGC STEMTCETFN DFVQLFKQVT SAEHLKDQTV YIVLDKAEYL RDMEANLLPG LLRLQELTDR NVTVIFLSEI IWEKFRPNTG CFEPFVLYFP DYSIGNLQKI LSHDHPPEYS ADFYAAYINI LLGVFYTVCR DLKELRHLAV LNFPKYCEPV VKGEAGERDT RKLWRNIEPH LKKAMQTVYL REISSSQWEK LQKDNTDPGQ LKGLSAYTHM ELPYYSKFIL IAAYLASYNP ARTDKRFFLK HHGKIKKTNF LKKHEKTSNH LLGPKPFPLD RLLAILYSIV DSRVAPTANI FSQITSLVTL QLLTIVGHED QLNGPKYKCT VSLDFIRAIA RMVNFDIKY LYDFL
Source	E.coli
Target Names	Orc5
Protein Names	Recommended name: Origin recognition complex subunit 5
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.