



# Recombinant Mouse Mitotic checkpoint protein BUB3 (Bub3)

<b>Product Code</b>	CSB-EP893001MO-B
<b>Abbreviation</b>	Bub3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9WVA3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTGSNEFKLN QPPEDGISSV KFSPNTSQFL LVSSWDTSVR LYDVPANSMR LKYQHTGAVL DCAFYDPHTA WSGGLDHQLK MHDLNTDQEN LVGTHDAPIR CVEYCPEVNV MVTGSWDQTV KLWDPRTPCN AGTFSQPEKV YTLVSVGDRL IVGTAGRRVL VWDLRNMGYV QQRRESSLKY QTRCIRAFPN KQGYVLSSIE GRVAVEYLDP SPEVQKKKYA FKCHRLKENN IEQIYPVNAI SFHNIHNTFA TGGSDGFVNI WDPFNKKRLC QFHRYPTSIA SLAFSNDGTT LAIASSYMYE MDDTEHPEDG IFIRQVTDAE TKPKST
<b>Source</b>	E.coli
<b>Target Names</b>	Bub3
<b>Protein Names</b>	Recommended name: Mitotic checkpoint protein BUB3 Alternative name(s): WD repeat type I transmembrane protein A72.5
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.