



Recombinant Putative T-box protein 38 (tbx-38)

Product Code	CSB-BP893180CXY
Abbreviation	tbx-38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XVD5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MYSCTSPSGI KVSLSTPEIW EEFYPKTEMI VTRNRGRVIF PHLDYIIKGL DPGSLYSIYI HLERVDGIKY KFDAGEWKEF AKGDPILPIQ YKEHPRGKRT GAEWMSEPVS FAHIKITNNP EIKDQKVILV QSMHKHIPVV TVKQVRHYKT GYQEDFSGEQ FRLEATEFMV VTAYQNEILK NLKVHHNKFA SGFRSNGKRR LSSDSENSEN SPPKRSKLV TPTISPQIDL PQQTPYYFNQ NFVAPQNYQP QFASAQNYNF EVQNAQLAW NMYQKQYEF WWQQQQMMMP GQPQELKNEF QSL
Source	Baculovirus
Target Names	tbx-38
Protein Names	Recommended name: Putative T-box protein 38
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.