



Recombinant Putative T-box protein 38 (tbx-38)

Product Code	CSB-MP893180CXY
Abbreviation	tbx-38
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9XVD5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MYSCTSPSGI KVSLSTPEIW EEFYPKTEMI VTRNRGRVIF PHLDYIIKGL DPGSLYSIYI HLERVDGIKY KFDAGEWKEF AKGDPIPIQ YKEHPRGKRT GAEWMSEPVV FAHIKITNNP EIKDQKVILV QSMHKHIPVV TVKQVRHYKT GYQEDFSGEQ FRLEATEFMV VTAYQNEILK NLKVHHNKFA SGFRSNGKRR LSSDSENSEN SPPKRSKLVV PPTISPQIDL PQQTPYYFNQ NFVAPQNYQP QFASAQNYNF EVQNNACLAW NMYQQKQYEF WWQQQQMMMP GQPQELKNEF QSL</p>
Source	Mammalian cell
Target Names	tbx-38
Protein Names	Recommended name: Putative T-box protein 38
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>