



Recombinant BTB/POZ domain-containing protein Y57A10B.3 (btb-14)

Product Code	CSB-MP893183CXY
Abbreviation	btb-14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XWH8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	TTELVEYTR GDIEFDSSNG EKIIDISTTK GIRCEWKVSS DYYTACFKWK FIWDQVDKYE VAGFSGQITV NYTTDSEGQK TRTVKVNLTG PGGEIWYNVS RSRYLSTYYA SYNITLEPQK RHMEAMELDK YFAPVDDRDA VLIVEGKKLH VSSCFLSFHS TYFHDLFEYD NSTSLLNIEI PVEGVSYEDL GLLLSVVHST ATFPNDGNSK KLELASQFQ TPYVLGLVEN HLLTNTFAWN ETLMLLADKY GLMRLLGKSI RRIDSIEKAR QLKDSVFLEL SDSAKVKVAY HMIQLL
Source	Mammalian cell
Target Names	btb-14
Protein Names	Recommended name: BTB/POZ domain-containing protein Y57A10B.3
Expression Region	22-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.