



Recombinant *Drosophila melanogaster* F-box/SPRY domain-containing protein 1 (Fsn)

Product Code	CSB-MP893237DLU
Abbreviation	Fsn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V6L9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVDPVAALCN YNVLEVIFSY LELDDLSHCS QVCKSWYHFL NDENSDVVRW HCLNKLPKES LKSDLLASVS TYKTKLRAYF HAWSPNDCSR NVYIKPNGFT LHRNPVAQST DAARGKIGFR HGRHTWEVIW EGPLGTVAVI GISTKEAALQ CHGYVALLGS DDQSWGWNLV ENHLLHNGDM QGSYPLLNNA PKYQVGERIR VILDCEENTL SFEKNYEFLG VAFRGLPDKK LYPTVSAVYG NTEVSMVYLG TPLDG
Source	Mammalian cell
Target Names	Fsn
Protein Names	Recommended name: F-box/SPRY domain-containing protein 1 Short name=DFsn
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.