



Recombinant Drosophila melanogaster Protein PTOV1 homolog (CG13609)

Product Code	CSB-MP893255DLU
Abbreviation	CG13609
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VCB1
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEPVPEEARG VGDDEGELDC LTHMIGALLL QKEPENPENG DSKETIWSGE LEWEDAQILD QPKTLHTVQC KICSMVKEGQ PEINTENWPN KLKTQLIPKK VLGKIGEQLF KDARMVVFRS SQGEVLNLLI TAMSSGFAGC IHFPSNPNCN IKALILYSR DHQALVGFIP NNEDSFSERL QEILQGAKRK PGVKTPKQPQ PEEPPPVEED AIINELLWTG SLNWSTQASL EEPSISHKLE CSVYIAIKNG DPGISAEDWP TDMPMVLMP5 IYLGQFAGAF IKDSKLIILR STPGEEHDSL ASSMSAGSCG CARFSSEVVC KVIMLLYSSP RNAFLGFIPR DQANFVKRLR EVLDEHRQKA RNKE
Source	Mammalian cell
Target Names	CG13609
Protein Names	Recommended name: Protein PTOV1 homolog
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.