



# Recombinant Drosophila melanogaster RuvB-like helicase 1 (pont)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP893279DLU  |
| <b>Abbreviation</b>      | pont   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q9VH07   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Drosophila melanogaster (Fruit fly)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MKIEEVKSTV RTQRIAAHSH VKGLGLDEVG AAVHSAAGLV GQKAAREAAG<br>IVVDLIKSKK MAGRALLLAG PPGTGKTAIA LAIAQELGNK VPFCPMVGSE<br>VFSNEIKKTE VLMENFRRSI GLRIRETKEV YEGETELTP VETENPMGGY<br>GKTISNVVIG LKTAKGTKQL KLDPSIFDAL QKEKVEVGDV IYIEANSGAV<br>KRQGRSDTFA TFDLETEEY VPLPKGDVHK KKEVIQDVTL HDLDVANARP<br>QGGQDVLSMM GQLMKPKKTE ITDKLRMEIN KVVNKYIDQG IAELVPGVLF<br>IDEIHMLDLE TFTYLHKSLE SPIAPIVIFA TNRGRCVIRG TTDIVSPHGI<br>PLDLLDRLLI IRTLLYSTAD MEQIIKLRAQ TEGQLQEENA FTRLSEIGTS<br>STLRYAVQLL TPAHQMCKVN GRNQISKDDI EDVHSLFLDA KRSSKHLSEK<br>NNKFML |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | pont   |
| <b>Protein Names</b>     | Recommended name: RuvB-like helicase 1 EC= 3.6.4.12 Alternative name(s):<br>Dpontin Short name= Dpon Pontin  |
| <b>Expression Region</b> | 1-456  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.