



Recombinant *Drosophila melanogaster* Vacuolar protein sorting-associated protein 72 homolog (YL-1)

Product Code	CSB-EP893294DLU
Abbreviation	YL-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VKM6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MAASRSRRNN AGNKIAHLLN EEEEDDFYKT SYGGFQEDEE DKEYEQKDEE EDVVDSDFSI DENDEPVSDQ EEAPEKRRKR GVVNTKAYKE TKPAVKKETK ATPALHKKRP GGGVTKRRPR PRFTVLDSCR KSIRTSTAIK TQATKIRLKE LDDARKRKKK KVRVEDYMPT QEELLEAKI TEEENTKSLE KFQKMELEKK KSRPTKRTFS GPTIRYHSLT MPAMRKPTRG ANPAVDSKDL AGKCERTFVT IENDFNDKVF QSLFRHKAPP KASNGICPIT RLPARYFDPI TQQPYYSIQA FKILREAYYM QLEQQGGGSE QPELAKWLEW RKLVKENRLK ASAAASKNGD N
Source	E.coli
Target Names	YL-1
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 72 homolog Alternative name(s): Protein YL-1
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.