



Recombinant *Drosophila melanogaster* WW domain-containing oxidoreductase (Wwox)

| | |
|--------------------------|--|
| Product Code | CSB-EP893301DLU |
| Abbreviation | Wwox |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9VLU5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila melanogaster</i> (Fruit fly) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MIALPDTDSE DELPPGWEER ATDDGTVCYV NQQGKTSQWT HPRTGRSKRI TGELPLGWK YDEQKGKRFM FLNKETQQR T NVDPRLAFV EEP TQNVAQV RQRFDSCSTA LQVLHGKDLH GRTALITGAN CGIGYETARS LAHHGCEIIF ACRNRSSAEA AIERIAQERP AARSRCRFAA LDLSSLRSVQ RFVVEIKQSV SHIDYLILNA GVFALPYTRT VDGLETTFQV SHLSHFYLT L QLET LFDYKT RIIVLSSESH RFANLPVENL AVHHLSPPE KYWSMMAYNN AKLCNVLFAQ ELAQRWKQRG ISVFSLHPGN MVSSDL SRNY WFYRLLFAIV RPFTKSLQQA AATSIYCATA NELTGLSGLY FNNCFFCEPS KLSKSAALQQ QLWK LSENLI AELVEQEQH |
| Source | E.coli |
| Target Names | Wwox |
| Protein Names | Recommended name: WW domain-containing oxidoreductase EC= 1.1.1.- |
| Expression Region | 1-409 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |