



Recombinant *Drosophila melanogaster* WW domain-containing oxidoreductase (Wwox)

Product Code	CSB-YP893301DLU
Abbreviation	Wwox
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VLU5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MIALPDTDSE DELPPGWEER ATDDGTVCYV NQQGKTSQWT HPRTGRSKRI TGELPLGW EK YYDEQ GKRFM FLNKETQ QRT NVDPR LAF AV EEPTQ NV AQV RQRF DSCSTA LQVLHGKDLH GRTALITGAN CGIGYETARS LAHHGCEIIF ACRNRSSAEA AIERIAQERP AARSRCRFAA LDLSSLRSVQ RFVVEIKQSV SHIDYLILNA GV FALPYTRT VDGL ETT FQV SHLSHFYLT L QLET LFDYKT RIIVLSSESH RFANLPVENL AVHHLSPPE KYWSMMAYNN AKLCNVLFAQ ELAQRWKQRG ISV FSLHPGN MVSSDLSRNY WFYRLLFAIV RPFTKSLQQA AATSIYCATA NELTGLSGLY FNNCFCEPS KLSKSAALQQ QLWK LSENLI AELVEQEQH
Source	Yeast
Target Names	Wwox
Protein Names	Recommended name: WW domain-containing oxidoreductase EC= 1.1.1.-
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.