



Recombinant *Drosophila melanogaster* Putative tRNA pseudouridine synthase Pus10 (Pus10)

Product Code	CSB-MP893311DLU
Abbreviation	Pus10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VPK7
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKNQDLVEYL RSCGVCEICQ LRYLKARGAE YRNINETLQK LDVKLNENVE NEELVISEGS QPSKRARLGV CSTCLGLFSK DFQNESSSI LASDFAKYDC QKIVLAICLP MSLQLRQLAM WFALRRRFGP SIDENNPPDV PIKEAVKLIL HPIICEKLAK VYDANGLMIN IDLTHSLEGE EVEKLVKLNK EAFPAKASHQ KRIEISRGLL EKQYQPSKVK AETYEKYFQI PCSNVEESIK LTSIDLQGPL ICVAGRYRKL SRELSHTPW I LNGQRLMEDS IEEIIIRHVG PHFTEKLGKI TFMSSGREDV DVRCLGKGRP FVLEIPNAKR SCLTKDEAFA MERAVDASGM VSIHHLQVVP REELTHIKTG EEQKKKFYRA LCALQEPVSV KILEQLQISA SFNIQQKTPI RVLHRRPLHT RPRTIYSVTA KVHRGNPKAL IIDIVTQAGT YIKELVHGEF GRTPPSLSSI IGKPMDIQAL DVVGIDLDWP HEVDNSKHKE
Source	Mammalian cell
Target Names	Pus10
Protein Names	Recommended name: Putative tRNA pseudouridine synthase Pus10 EC= 5.4.99.- Alternative name(s): Coiled-coil domain-containing protein 139 homolog tRNA pseudouridine 55 synthase Short name= Psi55 synthase tRNA pseudouridylate s
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.