



Recombinant *Drosophila melanogaster* Ubiquitin fusion degradation protein 1 homolog (Ufd1-like)

Product Code	CSB-MP893326DLU
Abbreviation	Ufd1-like
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VTF9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MFHFSGFNMM FPEGRNFHAN YKCFVSMLP GNERTDVEKG GKIIMPPSAL DTLTRLNVEY PMLFKLTNVK KSRSSHAGVL EFVADEGKCY LPHWMMENLL LGEGDILNIE SVSLPVATFS KFQPHSTDFL DITNPKAVLE NALRNFACTL RGDVIAIKYN KKVYELCVLE TKPGNAVSII ECDMNVEFEA PVGYKDHSET QASGSGQQGA AGTVGGEIAG ATNAILEEVV ETFKGGVRL DGKKKKESQL ETPVVKKVLA RGVDPDYDFQF GLIRFDRNIR PISDRSQEDD AVAGNADASD AESFHGTGFS MKKTRK
Source	Mammalian cell
Target Names	Ufd1-like
Protein Names	Recommended name: Ubiquitin fusion degradation protein 1 homolog Short name= UB fusion protein 1
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.